



Sterile Respirator Solutions

Protection for your worker. Assurance for your product.

Bullard

It's your life and you're worth it™

Worker Protection

- Full hoods offer 1000 Assigned Protection Factor (APF)
- Loose Fitting Facepieces offer 150% more protection than N95 masks.
- Patented Inflatable Neck Cuff
- Tychem QC construction

For use with Bullard's Industry Leading Powered Air Purifying Respirators (PAPRs)*

- Intrinsically Safe protection available
- Audible Low Voltage Alarms
- 33% more airflow than NIOSH minimum
- Up to ten hours of run time

Product Assurance

- Gamma Irradiated
- Pre-assembled
- Double bagged
- 10⁻⁶ Sterility Assurance Level

*The PAPRs mentioned in this document are not gamma irradiated and therefore no claim is made regarding sterility assurance for this component.



Ordering Information

| Part # | Description | | | |
|---------------|--|---|---|---|
| S-20TICRSA-IR | Sterile Assembly includes: 20TICN Double Bib Hood 20RT Ratchet Suspension PA1BT Breathing Tube Fully Assembled Packaged 5 per box. Minimum Quantity: 100 assemblies (20 boxes) |  |  |  |
| S-RT1SA-IR | Sterile Assembly includes: RT1 Double Bib Hood PA1BT Breathing Tube Fully Assembled Packaged 5 per box. Minimum Quantity: 100 assemblies (20 boxes) | |  |  |
| S-20FLSA-IR | Sterile Assembly includes: 20LFL (Large) Facepiece PA20LFBT Breathing Tube Fully Assembled Packaged 5 per box. Minimum Quantity: 100 assemblies (20 boxes) | |  |  |
| S-20LFMSA-IR | Sterile Assembly includes: 20LFM (Medium) Facepiece PA20LFBT Breathing Tube Fully Assembled Packaged 5 per box. Minimum Quantity: 100 assemblies (20 boxes) | | | |
| S-20TICHSA-IR | Sterile Assembly includes: 20TICH Double Bib Hood for Use with Hard Hats PA1BT Breathing Tube Fully Assembled Packaged 5 per box. Minimum Quantity: 100 assemblies (20 boxes) | |  |  |