

The **SilvIR facemasks** are made from extremely soft clean materials for greater operator comfort. They have a low Delta P rating for superior breathability and excellent bacteria filtration efficiency. They are validated sterile to a 10<sup>-6</sup> sterility assurance level in accordance with ANSI/AAMI/ISO guidelines

**DESCRIPTION:**

Product: Sterile cleanroom face mask  
Part Number: A-4FM7S  
Color: White  
Construction: Pleated with Ties

**MATERIALS:**

Outerfacing: 100% polyester  
Filter Media: White melt-blown polyester  
Inner Facing: 100% polyester  
Nosepiece: Flat aluminum wire

**Packaging:**

Individually packaged.  
25 per double bag, 12 packs per case.



**TEST DATA:**

B.F.E.: >95% post sterilization  
Delta P: 2.2  
Particulate: >0.5 um/minute = 98  
\* Test performed on a 9" mask (larger)

**STERILIZATION:**

Location: Steris, Libertyville, IL  
Method: Gamma Irradiation  
Validation: AAMI/ISO Method 1  
S.A.L.: 10<sup>-6</sup>  
Minimum Dose: 28.6 kGy  
Maximum Dose: 44.5 kGy