



NEW

Kimberly-Clark* Spring Green & Kimberly-Clark* Forest Green Nitrile Exam Gloves

Sustainable gloves – from start to finish

Looking for gloves that provide dexterity, strength, protection and a way to reduce your solid waste stream?


Try **NEW Kimberly-Clark* Spring Green and Kimberly-Clark* Forest Green Nitrile Gloves**. Comfortable as well as sustainable, the gloves:

- Reduce waste
- Save space
- Are recyclable

After use, these gloves won't go to waste. You can recycle them through **RightCycle* by Kimberly-Clark Professional***, the first large-scale recycling program for non-hazardous lab glove waste.

RightCycle* will turn your used gloves into the raw materials to create eco-friendly consumer products and durable goods.

It All Adds Up to a sustainable solution for a greener, safer, more productive lab.



All that Lab Waste Adds up Too

29%

29% of a university's solid waste comes from the Lab & Research Buildings.

22%

Lab gloves contribute 22% of the waste generated in a research lab.



Exceptional Workplaces*

Make your workplace exceptional

APEX is a unique employee training program that helps facilitate the changeover to the Kimtech® Brand of products. By effectively engaging employees, it ensures protection of the scientific process and enhances efficiency and performance.

NEW Kimberly-Clark* Spring Green and Kimberly-Clark* Forest Green Nitrile Exam Gloves

- Strong wet and dry grip to reduce potential for slips
- Textured fingertips for enhanced grip
- Beaded cuffs for added strength and donning ease
- Free of natural rubber latex – reducing the potential for glove-associated allergic reactions
- Superior quality – AQL exceeds the ASTM 2.5 standard for critical defects (pinholes)
- Food handling certified - Spring Green only
- Recyclable through RightCycle* by Kimberly-Clark Professional*

It All Adds Up: Our Gloves Increase Space, Reduce Waste and are Recyclable

1 RightCycle*

Our innovative recycling program that adds up to zero landfill and more efficient laboratories.



2 Space and Waste Reduction

Environmentally responsible packaging of our nitrile gloves that adds up to space and waste savings.



Environmentally Responsible Packaging (Compared with leading competitive gloves)

Kimberly-Clark* Spring Green Nitrile Exam Gloves – 2,000 count

100% More gloves per case
34% More storage space
15% Less waste (by weight)

AQL 0.65 exceeds ASTM 2.5 standard for critical defects (pinholes)

Kimberly-Clark* Forest Green Nitrile Exam Gloves – 2,000 count

100% More gloves per case
37% More storage space
51% Less waste (by weight)

AQL 1.5 exceeds ASTM 2.5 standard for critical defects (pinholes)

For more information on RightCycle* by Kimberly-Clark Professional*, visit www.kimtech.com/rightcycle

Kimberly-Clark* Spring Green and Kimberly-Clark* Forest Green Nitrile Exam Gloves

Description	Fingertip thickness (mils)	Color	XS	S	M	L	XL	XXL	Gloves/Boxes	Total/Case
Spring Green Nitrile Exam Gloves	4.7	Spring Green	43437	43438	43439	43440	43441	—	200/10 (XL=180/10)	2,000 (XL=1,800)
Forest Green Nitrile Exam Gloves	3.5	Forest Green	43443	43444	43445	43446	43447	43448	200/10 (XXL=180/10)	2,000 (XXL=1,800)

Our Guarantee

Your total satisfaction means everything to us. If, for any reason, our products do not meet your expectations, Kimberly-Clark will reimburse you* for your initial purchase, via FREE product, for up to \$1,000. For more information on KIMBERLY-CLARK PROFESSIONAL®, visit us online at www.kcprofessional.com.

* Guarantee extended to consuming end-user accounts only.



Exceptional Workplaces®

Kimberly-Clark warrants that its products (1) comply with K-C's standard specifications as of the delivery date to K-C's authorized distributors/direct purchasers (2) comply with all K-C labeling representations; and (3) are manufactured in compliance with all applicable federal, state, and local laws in effect at the time and place of manufacture of the products. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. K-C is not liable for any kind of special, incidental, or consequential damages. K-C's liability for breach of contract, tort or other cause of action shall not exceed the product purchase price. Purchasers and users are deemed to have accepted the above warranty and limitation of liability, and cannot change the terms by verbal agreement or by any writing not signed by K-C. To the extent required by applicable law, K-C does not limit its liability for death/injury resulting from K-C's negligence.