

Description

These pre-sterilized Suite Hose and Tubing Clamp are designed to be lightweight but used for high pressure applications. Its user-friendly operation allows for quick swing-away lockdown and can be added to existing assemblies without disturbing the process fluid path. The clamp is triple bagged and gamma irradiated for use in your most critical manufacturing environments.

Product Information

- Part Number:
 - V-BFS-101A-IR
- Packaging:
 - 1 clamp per double bag, 25 per pack, 4 packs per case



Performance

Test

Gamma Irradiation
Tensile Strength
Modulus of Elasticity
Flexural Strength

Result

25.0 – 50.0 kGy
180 – 200 MPa
9.50 – 10.5 GPa
270 – 290 MPa