

Sterile Polyethylene Trash Bags: Available in 3 stock sizes. Manufactured from high clarity polyethylene in an ISO Class 4 Cleanroom and validated sterile to a 10^{-6} sterility assurance level in accordance with ANSI/AAMI/ISO guidelines. Double bag packaged for aseptic entry into sterile core. Custom sizes available for other sterile bag needs. **MADE IN THE USA**

X-Stand: Autoclavable and chemical resistant 304 stainless steel X-stand holds 33x40 and 36x48 bags secure with visual assurance. **MADE IN THE USA**

Product:
 Sterile Polyethylene Trashbags

Part #:	Size:
PB2436-4S	24x36, 4ml
PB3340-4-IR	33x40, 4 ml
B-3648-1S	36x48, 1 ml

E-SS3340ST	Autoclavable X-stand
------------	-------------------------

Packaging:
 25 per double bag, 4 packs per case for the 4 mil. 10 packs per case for the 1 ml.



Testing:

Test	Method	Result
Sterility	ANSI/AAMI/ISO 11137	Validated Sterile to a 10^{-6} Sterility Assurance Level

- Also available: sterile polyethylene bucket liners!

Rev. 5-15