



VELTEK ASSOCIATES, INC.

ASEPTI-CLEANSE[®]

**Hands-Free DECON-AHOL[®] WFI
FORMULA/DECON-HAND[®] Dispenser**



**15 Lee Blvd.
Malvern, PA 19355
T: 610-644-8335
F: 610-655-8336
<http://www.sterile.com>**

Product Overview

VAI offers a hands-free dispensing system for DECON-AHOL® WFI Formula and DECON-HAND®. The dispensing system is the most advanced infrared sensor dispensing system available in the pharmaceutical and biotechnology industries. Developed to meet the requirements of cGMP clean room operations, ASEPTI-CLEANSE® is sealed to protect the internal electronics against cleaning and disinfecting. The unit mounts directly to glass or walls, making the dispenser an excellent choice for gowning rooms or aseptic manufacturing areas.

Quality and Manufacturing

- Is a sealed unit that protects the coated internal electronics
- Is photo-eye operated
- Just place your hand underneath and it dispenses a pre-measured dose to the hand without contact of the person to the unit
- Can be adjusted to dispense 1,3, or 5 mls
- Is designed in a dual power configuration of power supply; it can be powered by 4 D cell batteries or connected to a 110 volt receptacle in the same unit; In the electronically mounted installation, one has the backup of the battery power if a failure in the power line occurs due to abnormal situation
- Mounts directly on to glass or walls
- Is water resistant in design
- Is small; measuring 12 inches long x 5 inches wide x 2 inches high
- Has keyless opening and closing of the unit to change the contents

Order Number	Description	Qty/Cs
DEC-301	ASEPTI-CLEANSE Dispenser for 32oz Bottles	1
DECWFI-BOT-01	DECON-AHOL WFI Non-Sterile 70% IPA Bottles, 32oz Bottle	12
DECWFI-BOT-02	DECON-AHOL WFI Sterile 70% IPA Bottles, 32oz Bottle	12
DH-09	DECON-HAND Non-Sterile, 32oz bottle	12
DH-10	DECON-HAND Sterile 32oz bottle	12
DEC-301STAND	ASEPTI-CLEANSE Stand	1
DEC-301TRAY	ASEPTI-CLEANSE Drip Tray	1
DEC-301STANDSS	ASEPTI-CLEANSE Stand – Stainless Steel	1

