



KIMBERLY-CLARK* PURPLE NITRILE* Exam Gloves

Put the Science of Protection to Work in Your Lab

Your employees are your most valuable asset. That's why they deserve the worry-free protection of KIMBERLY-CLARK* PURPLE NITRILE* Exam Gloves. Our gloves guard against incidental splash exposure. The gloves are made of a highly durable nitrile polymer that eliminates the risks and costs associated with type 1 natural rubber latex allergies.

KIMBERLY-CLARK* PURPLE NITRILE* Exam Gloves are the gloves of choice in leading labs because they are comfortable to wear and provide excellent tactile sensitivity. When users see the distinctive color purple, they know they can concentrate on the task at hand, protected from known and unknown risks.



NEW SMARTPULL* Dispenser Box

New SMARTPULL dispensing incorporates 2 separate openings on the box. The first, smaller opening is used when the box is full and helps control dispensing to lessen waste from dropped gloves or multiple dispensing. When the box is half empty, the second, larger opening allows easier access to the gloves. In addition, the dispensing box also provides up to a 33% waste savings.*



Exceptional
Workplaces*

KIMBERLY-CLARK* PURPLE NITRILE* Exam Glove Features and Specifications

- Latex-free, powder free
- Beaded cuff
- Ambidextrous
- Textured fingertips
- AQL 1.5
- 9.5" & 12" lengths
- Excellent dexterity and tactile sensitivity
- Comfortable fit
- Cleared for use in chemotherapy
- PURPLE NITRILE* and PURPLE NITRILE-XTRA* Exam Gloves meet NFPA Standards 1999: 2008

Product Specifications			Physical Properties	
Gauge Thickness Measurements	MM	MIL	Before Aging	
Middle Finger	.15	5.9	Tensile Strength	21 MPa
Palm	.12	4.7	Ultimate Elongation	550%
Cuff	.09	3.5	After Aging	
Average Length	242mm (9.5")		Tensile Strength	21 MPa
	XTRA* 305mm (12.0")		Ultimate Elongation	500%
Quality Standards				
Exceeds current ASTM 6319 standard for critical defects (AQL 2.5). AQL for critical defects is 1.5.				
Manufactured in accordance with Quality System ISO 9001.				

Best Protection

When you're in the market for protection, turn to the proven leader. KIMBERLY-CLARK* PURPLE NITRILE* Exam Gloves provide peace of mind for your lab personnel.

KIMBERLY-CLARK* PURPLE NITRILE* Exam Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	Gloves/Boxes	Total/Case
KIMBERLY-CLARK* PURPLE NITRILE* Exam Gloves 9.5" Ambi	Purple	55080	55081	55082	55083	55084	100/10 (XL - 90/10)	1000 (XL - 900)
KIMBERLY-CLARK* PURPLE NITRILE-XTRA* Exam Gloves 12" Ambi	Purple	55090	50601	50602	50603	50604	50/10	500

Specific chemical resistance data can be found at: www.kc-safety.com/chemicalbarrierdata

Good for your business, good for the planet.

RIGHTCYCLE* from KIMBERLY-CLARK PROFESSIONAL* is an innovative program that helps you mitigate waste and cross-contamination issues in current processes and reach Corporate Social Responsibility (CSR) and Sustainability goals.

No more downcycling or upcycling. RIGHTCYCLE* makes it easy to recycle previously hard-to-recycle products like cleanroom garments and gloves into a variety of useful, eco-friendly products.

For more information, ask your distributor sales professional or contact KIMBERLY-CLARK PROFESSIONAL* directly at 800-255-6401.



Our Guarantee

Your total satisfaction means everything to us. If, for any reason, our products do not meet your expectations, Kimberly-Clark will reimburse you† for your initial purchase, via FREE product, for up to \$1,000. For more information on KIMBERLY-CLARK PROFESSIONAL*, visit us online at www.kcprofessional.com, contact your Kimberly-Clark Sales Representative, or call us at 1-888-346-GOKC (4652).

† Guarantee extended to consuming end-user accounts only.

Kimberly-Clark warrants that its products (1) comply with K-C's standard specifications as of the delivery date to K-C's authorized distributors/direct purchasers and are warranted for the following periods from end-user's date of purchase (verified by valid sales receipt) (a) five years for BALDER* Technology auto-darkening filters; (b) two years for all other auto-darkening filters; and (c) one year for powered air-purifying respirators; (2) comply with all K-C labeling representations; and (3) are manufactured in compliance with all applicable federal, state, and local laws in effect at the time and place of manufacture of the products. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. K-C is not liable for any kind of special, incidental, or consequential damages. K-C's liability for breach of contract, tort or other cause of action shall not exceed the product purchase price. Purchasers and users are deemed to have accepted the above warranty and limitation of liability, and cannot change the terms by verbal agreement or by any writing not signed by K-C. To the extent required by applicable law, K-C does not limit its liability for death/injury resulting from K-C's negligence.