

# MICROFLEX®

## ▲ LONG CUFF NITRILE EXAM GLOVES



TEAL  
NITRILE  
NON-STERILE  
AMBIDEXTROUS  
FULLY TEXTURED  
BEADED CUFF  
FOR SINGLE USE ONLY  
NON-LATEX

### POWDER-FREE EXAMINATION GLOVES

The long cuff and increased thickness provides an additional level of hand protection for the user. Nitrile offers superior puncture resistance. A textured finish allows for a good grip.

- Longer length and thicker for added protection



PRODUCT CODE		PACKAGING	
S	N891	50	Gloves/Box
M	N892	10	Boxes/Case
L	N893		
XL	N894		

GLOVE PROFILE		ASTM D 6319	N89
Length	<i>mm (inches)</i>	min. 220 (8.7)	300 (11.8)
Cuff Thickness	<i>mm (mils)</i>	n/a	0.10 (3.9)
Palm Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.15 (5.9)
Finger Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.18 (7.1)
Tensile Strength	<i>Before Aging</i>	min. 14	27
	<i>After Aging</i>	min. 14	14
Elasticity/Elongation	<i>Before Aging</i>	min. 500	500
	<i>After Aging</i>	min. 400	400
Pinhole Defect Rate*	<i>(AQL)</i>	max. 2.5	1.5

\* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED  
INDUSTRIES