



VELTEK ASSOCIATES, INC.

DECON-Clean[®] WIPE

**Cleaning Agent/Residue Remover
Saturated Wipe**



**15 Lee Blvd.
Malvern, PA 19355
T: 610-644-8335
F: 610-655-8336
<http://www.sterile.com>**

Product Overview

DECON-Clean wipe RTU is a superior saturated wipe designed to remove residue left behind from disinfecting agents. It is used on both aseptic and non-aseptic wipe downs of filling and packaging machinery, stainless steel, lexan, polycarbonate, glass, and any critical surfaces that require cleaning. The unique and patented Asepti-Fill system and gamma irradiation processing assures sterility of the wipes is contained.

Use DECON-Clean wipes for areas that require the assurance of a surface that is free from cleaning agents and product residues.

Features and Benefits

- Available in both sterile and non-sterile 12" x 12" versions
- Sterile version, gamma irradiated 10⁻⁶ SAL
- Delivered with the lot specific Certificate of Analysis, Certificate of Irradiation, and Certificate of Sterility
- Sterility tested according to current USP compendium
- Cut using FocusEdge cutting technology
- Extremely low in particulate shedding and soluble extractable
- Laundered in ultra-filtered water of 0.45µ
- Dried in a ULPA filtered Class 1 air stream
- Packaged Quadruple Bagged using the ABCD Clean Room Introduction System™
- Packaged 20 wipes per pack that are folded for easy removal

Order Number	Description	Qty/Cs
VEL14-12X12-S-4002	DECON-CLEAN Wipe 12" x 12" Sterile	200
VEL14-12X12-4003	DECON-CLEAN Wipe 12" x 12" Non-Sterile	200