



## SOFTWEAR<sup>®</sup> Nitrile Exam Gloves

Powder-Free  
9.5" Cuff  
4.5 mil (palm)  
Medical/Exam Grade  
White Color

*The ultimate in fit, feel and comfort.*

A soft and stretchy white nitrile exam glove that offers excellent durability and tactile sensitivity. This glove is fully textured for enhanced gripping power.

Cat#	Size
N791	Small
N792	Medium
N793	Large
N794	X-Large

Packaging: 100 Gloves/Box, 10 Boxes/Case

### SPECIFICATIONS

- Finish: Textured
- Type: Ambidextrous / Non-Sterile
- Cuff: Beaded
- Thickness: 4.5 mil (palm); 4.0 mil (cuff)  
5.5 mil (finger);  $\pm 0.01$ mm
- Tensile Strength: 20.5 Mpa (min)
- Elasticity/Elongation: 700% (min)
- 1.5 AQL
- Complies with FDA Biocompatibility Guidance for Medical Devices
- 510K/Approved for medical use
- ISO 9001 Certified QMS