

Suite[®] STERILE RESPIRATOR SOLUTIONS

Protection for your worker.
Assurance for your product.



Superior Worker Protection

- Loose-fitting facepieces provide 150% more protection than N95 masks
- Patented air delivery systems help ensure positive pressure
- DuPont™ Tychem® QC construction offers superior protection



Key Features

- Built-in soft suspension
- 25 APF (Assigned Protection Factor)
- No fit testing required
- Wide angle lens for increased visibility
- Filters 99.97% of particulates from exhaled air
- Chemical, water, and moisture resistant
- Constructed of DuPont™ Tychem® 2000
- NEW HEPA exhalation panel



Sterile Product Confidence

- Gamma irradiated
- Pre-assembled and Double bagged
- Validated sterile



Ordering Information

Part Number	Description	Packaging
• S-20LFMHE-S	Sterile 20LF Facepiece, Standard Width, Medium, HEPA Exhalation Panel, Breathing Tube Attached	5 per box
• S-20LFLHE-S	Sterile 20LF Facepiece, Standard Width, Large, HEPA Exhalation Panel, Breathing Tube Attached	5 per box
• S-20LFXLHE-S	Sterile 20LF Facepiece, Standard Width, X-Large, HEPA Exhalation Panel, Breathing Tube Attached	5 per box