

**Description**


Suite<sup>®</sup> wipers are engineered for high-performance. Uniquely manufactured from synthetic high-density ultra fibers, Suite wipers are non-particulating, snag resistant, and offer unrivaled absorbency. They deliver maximum performance and durability when used as part of your cleaning and/or disinfecting protocols.

**Featuring:**

Double welded seam to reduce particle generation in comparison to traditional knife cut or single welded seam.

**Product Information**

- Part Numbers:
  - Non-Sterile: W-SC12-25
  - Sterile: W-SC12-5S
- Size: 12" x 12"
- Color: White
- Packaging:
  - Non-Sterile: 200 wipers/case (25 wipers/bag, 8 bags/case)
  - Sterile: 200 wipers/case (5 wipers/bag, 4 bags/pack, 10 packs/case)



**Product Performance**

Test	Method	Result
Particles	IEST-RP-CC03	Acceptable for ISO Class 5 use
Sterility	AAMI TIR33:2005	Validated Sterile to a 10 <sup>-6</sup> Sterility Assurance Level

**Reduction and Transfer:** Wiping with Suite material has been demonstrated to \*effectively remove surface contamination. In addition, the carry-over (cross contamination potential) after wiping with the Suite material is minimized as contaminants are fully retained by the material.

\*Lab testing available upon request