




Polyester Wipers 9"x9"

Specification Sheet

Description

The BenchmarkGMP Polyester Wipers are made from a 100% continuous-filament, double-knit polyester fabric. They are cleanroom packed and laundered with product traceability. Their ultrasonic sealed edge means they are extremely low particle and fiber generating. They are an ideal economic and quality wiper choice.



Applications

- Ideal for ISO 4 cleanrooms and above
- Cleaning and disinfecting in all areas of cleanroom and production

Specifications

Type	Method	Results
Fibers	Orbital Shake Test IEST-4.3	< 400 counts/m ²
LPC	Orbital Shake Test IEST-4.3	6.0 x 10 ⁶ /m ²
Organic Contaminates		Does not contain Silicon, Amides or DOP

Part Number	W-PKUW-99		
Size	9"x 9"		
Packaging	150 wipers per bag, 10 bags per case		
Material	100% Virgin Grade Polyester (75 Denier)		
Edge Model	Ultrasonic Sealed Edge		
Basis Weight	135 g/m ²		
Sorptive Capacity	380 ml/m ²		
Sorptive Rate	< 1 second		
Ions Contamination (ppm)	Calcium: 0.09	Sodium: 0.20	Potassium: 0.03
NVR (g/m ²)	DI-Water: ≤ 0.02	IPA: ≤ 0.03	

