

ROUGING REMOVAL

SUITE® SPONGES, WIPERS & CITRISURF



DESCRIPTION

Current disinfection procedures leave disinfectant residues on cleanroom surfaces which lead to facility appearance issues and corrosion stains on stainless steel equipment.

Suite® products are high performance engineered cleaning products made from unique synthetic high-density ultra-fibers which are non-particulating, snag resistant, offer unrivaled absorbency, and thoroughly clean hard surfaces in any environment. The CitriSurf® citric acid cleaning and passivation process is based on cutting edge technology using safe and environmentally friendly citric acid to remove the free iron from the surface of stainless steel.

Combining the power of Suite Sponge with Citrisurf Solutions allows you to remove even the toughest stainless steel rouging.

PRODUCT INFORMATION

Suite Sponge Part Numbers:

1. Sterile: Y-GL-SP47-1S & Y-GL-SP47-10S
2. Non Sterile: Y-GL-SP47-10

Suite Wiper Part Numbers:

1. Sterile: W-TTC1212-CS & W-TTC1212-25CS
2. Non Sterile: W-TTC1212-C

Citrisurf Part Numbers:

77 Plus: Liquid Spray for flat surfaces

1. C-CMP12-0077P. 12 x 22 oz. Spray Bottles/Case

2310: Gel Spray for vertical, horizontal, and uneven surfaces

2. C-CMP12-2310. 12 x 22 oz. Spray Bottles/Case



BP-STE-0008-A

