

CLEANROOM WIPERS

SUITE® SC6 WIPERS



DESCRIPTION

Suite® wipers are engineered for high-performance. Uniquely manufactured from synthetic high-density ultra fibers, Suite® wipers are non-particulating, snag resistant, and offer unrivaled absorbency. They deliver maximum performance and durability when used as part of your cleaning and/or disinfecting protocols.

- Double welded seam to reduce particle generation in comparison to traditional knifecut or single welded seam
- The wiper's textured pattern provides improved surface cleaning capabilities.



PRODUCT INFORMATION

Part Numbers:

- **Non-Sterile:** W-SC6-5
- **Sterile:** W-SC6-5S

Size: 12" x 12"

Color: White

Packaging:

- **Non-Sterile:** 250 wipers/case (5 wipers/bag, 50 bags/case)
- **Sterile:** 250 wipers/case (5 wipers/bag, 5 bags/pack, 10 packs/case)



PRODUCT PERFORMANCE

Type	Method	Result
Particles	IEST-RP-CC03	Acceptable for ISO Class 5 Use
Sterility	AAMI TIR33:2005	Validated Sterile to a 10 ⁻⁶ Sterility Assurance Level

Reduction and Transfer: Wiping with Suite® material has been demonstrated to *effectively remove surface contamination. In addition, the carry-over (cross contamination potential) after wiping with the Suite® material is minimized as contaminants are fully retained by the material.

*Lab testing available upon request.

BP-STE-0008-A

