

SUITE® SINGLE USE BETA BAG



DESCRIPTION

The Suite Single-Use Beta Bag is a safe, efficient, and cost-effective solution for aseptic production. This single-use system provides contamination-free transfer of sterile components into isolation or removal of waste components out of isolation. Used in conjunction with Alpha port systems, sealed containment and sterility of the environment are maintained.

The Single-Use Beta Bag system consists of Tyvek® and high-density polyethylene (HDPE) bag layers connected to a Beta flange which mates the entire assembly to the Alpha port on an isolator. Once filled with components, these bags can be autoclaved and connected for filling line production. Bag customization is also available to meet specific application needs based on component, processing, and production parameters.

In addition to its proprietary design the Suite Single-Use Beta Bag flange is compatible with other Alpha port designs currently available. This allows manufacturers more flexibility with additional beta bag options available at their disposal.



FEATURES & BENEFITS:

- Fits 190 mm diameter Alpha Port
- Available in both autoclavable and gamma-compatible versions
- Available gamma irradiated for waste removal
- Custom gamma irradiated kitting solutions available
- Removable, reusable handle available
- Allows for increased flexibility in production
- Single-use alternative to traditional stainless steel transfer containers maximizes productivity by reducing downtime needed for cleaning
- Tyvek® and HDPE layers form sterile barrier for autoclave sterilization process



BP-STE-0001-B



SINGLE USE BETA BAG



COMPONENTS



SINGLE-USE BETA BAG

BETA COVER WRENCH



Use the Beta Wrench to open/close the Beta Bag when placing contents inside

REMOVABLE HANDLE



Use the Removable Handle when docking/undocking the Beta Bag to the Rapid Transfer Port

Product Number	Description	Port Size	Bag Size	Quantity
B-SBB-50547	Tyvek/HDPE for Autoclave	190	25L	2/Case
B-SBB-50990-8IR	Tyvek/HDPE Gamma-Irradiated	190	25L	2/Case
B-SBB-50990	Tyvek/HDPE for Gamma	190	25L	2/Case
S-49500	Removeable Handle System	190	-	1/Case
S-41115	Beta Cover Wrench	190	-	1/Case

