

# CLEANROOM EQUIPMENT

## SUITE® SUDS FRAME



### DESCRIPTION

The Suite® SUDS frame utilizes the dynamic fiber technology of the Suite Touchless mops while being extremely light-weight and durable. This gamma irradiated disposable mop frame eliminates the time and cost of cleaning and autoclaving equipment and is also comprised of 100% recyclable polypropylene for sustainability. The frame hinge features a double ratcheting system which allows the operator to maintain the mop in the position needed to achieve their ideal mopping configuration.

### PRODUCT INFORMATION



#### Part Numbers:

- **Sterile:** Y-SUDS-2-20S
- **Non-Sterile:** Y-SUDS-2-20
- **Handle 37" - 71" Extension:** Y-SUITE-HDL
- **Touchless Mop Pad:**
  1. Y-TL-517-25
  2. Y-TL-517-25S
- **Touchless Floor Mop:**
  1. Y-TL-FM-25
  2. Y-TL-FM-25IR

#### Material Composition:

- 100% recyclable polypropylene

**Packaging:** 1 frame per package



100% Recyclable



BP-STE-0001-B

