

# CLEANROOM SPONGES

## STERILE SUITE® SPONGES



### DESCRIPTION

Suite® products are high performance engineered cleaning products made from unique synthetic high-density ultra-fibers which are non-particulating, snag resistant, offer unrivaled absorbency and thoroughly clean hard surfaces in any environment achieving supreme performance and durability.

This validated sterile cleanroom hand sponge is comprised of polyurethane foam wrapped with a bound-seam sealed and hemmed cover made of the Suite® Ultra-Fiber textile technology.



### PRODUCT INFORMATION

#### Part Numbers:

1. Sterile: Y-GL-SP47-1S
2. Sterile: Y-GL-SP47-10S Non-Sterile: Y-GL-SP47

#### Packaging:

1. 1 sponge per bag, 10 bags per pack, 10 packs per case
2. 10 sponges per pack, 10 packs per case

**Size:** 4"x7"

**Weight:** 0.6 oz.

**Color:** White

**Ability to Pick-Up and Retain:** < 0.2 micron particles



### PRODUCT PERFORMANCE

Type	Method	Result
Sterility	AAMI TIR33:2005	Validated Sterile to a 10 <sup>-6</sup> Sterility Assurance Level
Particles	IEST-RP-CC03	Acceptable for ISO Class 5 Use

**Reduction and Transfer:** Wiping with Suite® material has been demonstrated to \*effectively remove surface contamination. In addition, the carry-over (cross contamination potential) after wiping with the Suite® material is minimized as contaminants are fully retained by the material.

\*Lab testing available upon request.

BP-STE-0008-A

