

CLEANROOM EQUIPMENT

SUITE[®] SPRAYER



DESCRIPTION

The Suite Sprayer is a versatile and efficient solution for effective application of disinfectants, cleaning solutions, and sterilants. This gamma irradiated sprayer is designed with exceptional features to enhance usability and deliver consistent, high-quality performance.

FEATURES

- Adjustable nozzle allows for precise flow rate control (0.4-0.5 gpm)
- 48-ounce tank capacity for extended use
- Lightweight construction (0.75lb) reduces operator fatigue
- Double bagged within two case liners for Cleanroom Grade A and B area use
- Pump-driven system ensures constant, even pressure for uniform application
- Ergonomically designed trigger spray for comfortable handling
- Versatile spray patterns adapt to different surfaces and areas
- Vertical stream of 23 ft, horizontal 25 ft



Part Number

- S-RTS-1098-CIR

Packaging

- 15 bottles per case

Sterility

- Gamma irradiation

BP-STE-006-C

