

2-PLY POLYESTER WIPERS 12"X12"

SPECIFICATION SHEET



DESCRIPTION

The BenchmarkGMP 2-Ply Polyester Wipers are made from a 100% continuous-filament, double-knit polyester fabric. They have a clean and durable ultrasonic bonded 2-layer construction which gives them higher absorbency compared to 2 single wipers. Their border sealed edge means they are the least particle and fiber generating edge compared to other edge sealings.

APPLICATIONS

- Ideal for ISO 4 cleanrooms and above
- Cleaning and disinfecting in all areas of cleanroom and production



SPECIFICATIONS

Type

Fibers
LPC
Organic Contaminates

Method

Orbital Shake Test IEST-4.3
Orbital Shake Test IEST-4.3

Results

< 600 counts/m²
7.0 x 10⁶ particle/m²
Does not contain Silicon, Amides or DOP

Part Number	W-Q2PLYBS-1212 W-Q2PLYBS-1212-S
Size	12"x 12" 12"x 12" (Sterile)
Packaging	100 wipers per bag, 5 bags per case 20 wipers per bag, 5 bags per pack, 5 packs per case
Material	100% Virgin Grade Polyester (75 Denier)
Edge Model	Ultrasonic Border Sealed Edge
Basis Weight	220 g/m ²
Sorptive Capacity	680 ml/m ²
Sorptive Rate	< 1 second
Ions Contamination (ppm)	Calcium: 0.10 Sodium: 0.25 Potassium: 0.05
NVR (g/m²)	DI-Water: ≤ 0.02 IPA: ≤ 0.03