

CLEANROOM SOCKS

SPECIFICATION SHEET



DESCRIPTION

BenchmarkGMP Cleanroom Socks - the ideal solution to keep manufacturing facilities Annex 1 compliant. These cleanroom laundered and packaged socks are made of continuous filament polyester and nylon yarn. They minimize bioburden levels and reduce particulate contamination in critical environments. The mid-calf length and elastic cuff provide a comfortable and secure fit, while the polyester/nylon material ensures durability. Designed to fit all foot sizes, these lint-free socks prevent the risk of contaminating your cleanroom with dirty street socks.

Compliant with EU GMP Annex 1 regulations, the BenchmarkGMP Cleanroom Socks meet the requirements for entering changing rooms leading directly to Grade B and Grade C cleanrooms. They provide maximum coverage with facility suits, covering the full length of the arms and legs, and socks covering the feet, ensuring no outdoor clothing is brought into the changing rooms. Additionally, the socks reduce the risk of contamination in the gowning areas or processes, maintaining an aseptic environment.



SPECIFICATIONS

Type
Particles

Method
Helmke Drum Test
IEST-RP-CC003.4

Results
Category I

Part Numbers	A-BP-SOCK2-S/M, A-BP-SOCK2-L/XL, A-BP-SOCK2-2X
Part Numbers (Gamma Irradiated)	A-SOCK-S/M-IR, A-SOCK-L/XL-IR, A-SOCK-2X-IR
Size	Small-Medium, Large-XLarge, 2X-Large
Material	Polyester/Nylon
Color	White
Packaging	10 pairs per pack, 50 packs per case

V.7.BC

