

BENCHMARK GMP STERILE 9" FACEMASK WITH TIES

BENCHMARK GMP

PRODUCT INFO

Sterile
GMP Grade A compatible
ISO Class 4 & 5 compatible
9 inches length with malleable nose
strip

MATERIALS

TIE-TAPES

Non-woven Polypropylene

OUTER

Non-woven Spunbond Polyester

FILTER

Non-woven Meltblown Polyester

INNER

Non-woven Polypropylene/
Polyethylene

QUALITY ASSURANCE

Manufactured in a facility operating
under an ISO 9001:2015 quality
management system

Processed in a NEBB certified ISO Class
5 Cleanroom

STERILIZATION

Gamma irradiation, Minimum 25 kGy

SHELF LIFE

Three (3) years from date of
manufacture

STORAGE

Store in a dry, cool place (<40°C) away
from direct sunlight. Do not expose
open cartons to prolonged direct
fluorescent light

PACKING

One piece per sealed inner PE bag
50 inner bags per sealed outer PE bag
10 outer bags per lined carton (500
pieces)

SIZE

Universal



PHYSICAL PROPERTIES

LENGTH	WIDTH	COLOUR	STYLE
9" (+/-0.2")	3.7" (+/-0.2")	White	Tie-On

FILTRATION PROPERTIES

PARTICLE FILTRATION EFFICIENCY (PFE) @ 0.1 MICRON	94%
BACTERIAL FILTRATION EFFICIENCY (BFE) @ 0.3 MICRON	97%
DIFFERENTIAL PRESSURE (ΔP), MM H ₂ O/CM ²	6.8

PART NO

A-BP9FM4T-S