

SINCE 1925

KEYSTONE

SAFETY ✦ CLEANROOM PRODUCTS

BLUE BHD MUFFS



Used as an autoclavable cover, these Blue Heavy Duty (BHD) Muffs eliminate the need for time consuming wrapping with sheets of Bioshield™ material and autoclave tape.

TYPICAL APPLICATION & CHARACTERISTICS:

- Manufactured and packaged in ISO certified Cleanroom
- Autoclavable up to 135°C
- Conveniently packaged: 100 per heat-sealed, labeled poly bag
- Color: Blue
- Certificates of Conformance sent with all shipments
- Made in USA

MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	BHD MATERIAL
Basis Weight	Gm/m2	57
Dry Tensile Strength (MD)	grams	8,100
Dry Tensile Strength (CD)	grams	7,100
Elongation Dry	%	27.3
Air Permeability	liters / 100 sq. cm	110
Suter Water Head (pressure test)	mm	400
Dry Mullen Burst	Grams / sq.cm	2,150
Flammability	seconds	8

SINCE 1925

KEYSTONE

SAFETY ✦ CLEANROOM PRODUCTS

BLUE BHD MUFFS



KEYSTONE PART NUMBER	DIMENSIONS
MFNWB03	3" used for openings 1/2" to 1"
MFNWB04	4" used for openings 1" to 1 3/4"
MFNWB05	5" used for openings 1 3/4" to 2 1/2"
MFNWB07	7" used for openings 2 1/2" to 4"
MFNWB08	8" used for openings 3 1/2" to 4 3/4"
MFNWB09	9" used for openings 4" to 5 1/2"
MFNWB12	12" used for openings 5 1/2" to 8"
MFNWB16	16" used for openings 8" to 12"
MFNWB18	18" used for openings 10" to 14"