

CLASS 100 CONTROLLED ENVIRONMENT LATEX GLOVES

PFXT SERIES

CTI's PFXT Series are Class 100 Gloves. Like all our gloves, they are manufactured for outstanding strength, stretch and fit. They offer a textured surface and are double chlorinated reducing contact sensitizers and protein content.

PFXT Series gloves are ISO 5 Cleanroom compatible or higher for controlled environments. We offer them in both 10" and 12" long. They are ambidextrous, and have been tested to meet and/or exceed all standards and requirements. Available in sizes Small to X-Large.



PRODUCT FEATURES

- Textured Surface
- Low Protein and Contact Sensitizers
- Class 100 (ISO 5) Cleanroom Compatible or Higher for Environments
- Double Chlorinated
- Double Bagged

- Total Protein Content <50ug/g of Glove (ppm) Using Modified Lowry Method
- Leak Tested
- Meets ASTM D3578-05 standards
- Sampled in accordance with MIL STD-105E
- Complies with U.S. Government 21 CFR 177.2600 for food handling
- Qualified for Use in the Pharmaceutical Industry

SPECIFICATIONS & PACKAGING

Length: **10" and 12"** (+/-0.5)

Thickness: **6.5 mil** (+/-1)

Tensile Strength: **3400 psi/23 mpa**

Elongation: **745%**

Packaging: **Poly Bagged, 100 Gloves/Pack, 1000 Gloves/Case**

Cat No./Sizes:

10" Sizes:

PFXT101 – Small

PFXT102 – Medium

PFXT103 – Large

PFXT104 – X-Large

12" Sizes:

PFXT201 – Small

PFXT202 – Medium

PFXT203 – Large

PFXT204 – X-Large