

BENCHMARK GMP STERILE CLEANROOM SLEEVE COVERS

PRODUCT INFO

ISO Class 4 (Class 10) & above
GMP-Kleen material
Extra long 22 inches with thumb-loop
Sterilization via gamma irradiation
SAL 10⁻⁶
Latex-free with bound seams and
enclosed elastic

QUALITY ASSURANCE

Manufactured in a facility operating
under ISO 9001 certification

Processed in a NEBB certified
ISO Class 5 Cleanroom

SHELF LIFE

Three (3) years from date of
manufacture

STORAGE

Store in a dry, cool place (<40°C) away
from direct sunlight. Do not expose
open cartons to prolonged direct
fluorescent light

PACKING

1 pair per sealed bag,
10 pairs per outer bag,
10 outer bags per lined carton (100
pairs)

SIZING

UNIVERSAL

PART NO

A-BP-SLV20-100-S

COUNTRY OF ORIGIN: MALAYSIA



MATERIAL PERFORMANCE PARAMETERS

DESCRIPTION	UNITS	METHOD	SPECIFICATION	TYPICAL
Basic Weight	g/m ²	ISO 9073-1-1989	±2	60
Machine Direction strength	N/5cm		≥50	69
Machine Direction Elongation	%		≥100	109
Cross Direction Strength	N/5cm		≥25	43
Cross Direction Elongation	%		≥100	111
Hydrostatic Head	mmH ₂ O	ISO 811-1981	≥100	213
Air Permeability	l/m ² /S	GB/T5453	≥	2000
Surface Resistivity	Ω	ISO 9073-1-1989	/	10 ⁻⁷

CLEANLINESS LEVELS	UNITS	METHOD	SPECIFICATION	TYPICAL
Particle Shedding	≥ 0.5µm (counts/min)	IEST-RP-CC003.2	2000	1000