


PUREGRIP® Glass Bottles

www.foxxlifesciences.com



- Available in clear and amber options
- · Sizes range from 100 mL 20 L
- Highest Polypropylene cap temperature rating in the industry, 142°C
- Industry standard GL45 thread, compatible with Autofil™ Bottle-Top Filters and EZBio® Single-Use Cap Assemblies
- Large, visible writing surface on bottle for additional labeling; white graduations for identification of approximate volume; retrace code for identification
- BOROSIL® 3.3 Borosilicate Glass is mechanically strong, chemically resistant, and withstands thermal shock
- Bottle, VersaCap®, and Pouring Ring are autoclavable and complies with USP-VI

PUREGRIP® Round Glass Media Storage Bottles combine the leading quality of BOROSIL® 3.3 Borosilicate Glass with the patented advantages of the new Foxx GL45 VersaCap® to produce a true leap forward in laboratory glassware. Ideal for mixing, storage, and transportation of media, buffers, chemicals, and solvents, **PURE**GRIP® Bottles offer versatility across a wide variety of laboratory needs. And with Foxx's revolutionary GL45-threaded polypropylene VersaCap®, **PURE**GRIP® Bottles are substantially easier to hold securely, improving lab safety and user experience.

PUREGRIP® Bottles are autoclavable, reusable, and graduated in mL. A large writing space allows for visible labeling, and an extra frosted surface on the VersaCap® offers a convenient second location for writing. Manufactured from top-quality glass, **PURE**GRIP® Bottles are highly resistant to thermal shock. And with an industry-leading maximum temperature rating of 142°C, the VersaCap® stands up to even the toughest conditions.



2 Patented SECUREGRASP technology built into cap for ultimate ease of handling, especially with gloved hand

3 USP Class VI, Animal-Free Polypropylene VersaCap® and Pouring Ring for leakproof, drip-free experience

- 4 Sloped shoulders for smooth, natural pouring motion
- 6 Glass bead indicates full capacity level
- 6 Large, visible writing surface on bottle for additional labeling; white graduations for identification of approximate volume; retrace code for identification
- 7 Easy to read graduation marks



PUREGRIP® is a registered trademark of Foxx Life Sciences.





PUREGRIP® Glass Bottles

www.foxxlifesciences.com



	Product Number	Description	Total Capacity (mL)	Approx. Dimensions (Dia. x H)	Thread Type	Quantity Per Case
Clear	1501016-FLS	PUREGRIP® Reusable Glass Media Storage Bottles with VersaCap®	100	56 x 103	GL45	10
	1501021-FLS	PUREGRIP® Reusable Glass Media Storage Bottles with VersaCap®	250	70 x 141	GL45	10
	1501024-FLS	PUREGRIP® Reusable Glass Media Storage Bottles with VersaCap®	500	86 x 179	GL45	10
	1501029-FLS	PUREGRIP® Reusable Glass Media Storage Bottles with VersaCap®	1,000	101 x 229	GL45	10
	1501030-FLS	PUREGRIP® Reusable Glass Media Storage Bottles with VersaCap®	2,000	136 x 265	GL45	10
	1501031-FLS	PUREGRIP® Reusable Glass Media Storage Bottles with VersaCap®	3,000	158 x 298	GL45	2
	1501033-FLS	PUREGRIP® Reusable Glass Media Storage Bottles with VersaCap®	5,000	183 x 326	GL45	1
	1501038-FLS	PUREGRIP® Reusable Glass Media Storage Bottles with VersaCap®	10,000	235 x 408	GL45	1
	1501040-FLS	PUREGRIP® Reusable Glass Media Storage Bottles with VersaCap®	20,000	296 x 503	GL45	1
Amber	1519016-FLS	PUREGRIP® Reusable Amber Glass Media Storage Bottles with VersaCap®	100	56 x 103	GL45	10
	1519021-FLS	PUREGRIP® Reusable Amber Glass Media Storage Bottles with VersaCap®	250	70 x 141	GL45	10
	1519024-FLS	PUREGRIP® Reusable Amber Glass Media Storage Bottles with VersaCap®	500	86 x 179	GL45	10
	1519029-FLS	PUREGRIP® Reusable Amber Glass Media Storage Bottles with VersaCap®	1,000	101 x 229	GL45	10
	1519030-FLS	PUREGRIP® Reusable Amber Glass Media Storage Bottles with VersaCap®	2,000	136 x 265	GL45	10
	1519031-FLS	PUREGRIP® Reusable Amber Glass Media Storage Bottles with VersaCap®	3,000	158 x 298	GL45	2
	1519033-FLS	PUREGRIP® Reusable Amber Glass Media Storage Bottles with VersaCap®	5,000	183 x 326	GL45	1
	1519038-FLS	PUREGRIP® Reusable Amber Glass Media Storage Bottles with VersaCap®	10,000	235 x 408	GL45	1
	1519040-FLS	PUREGRIP® Reusable Amber Glass Media Storage Bottles with VersaCap®	20,000	296 x 503	GL45	1

*EZBio® Cap Assemblies are compatible with the GL45 VersaCap® and can be purchased separately



Rev. 08/19/2022 PN: DS00013-01

