

PRE-SATURATED 12" X 12" STERILE POLYESTER KNIT WIPES 70% IPA/30% WFI



SPECIFICATION SHEET

DESCRIPTION

These sterile pre-saturated double-knit polyester wipes are made from 100% pure, continuous filament, polyester. The durable construction makes it the ideal choice for cleanrooms and critical manufacturing environments. Cleanroom laundered, the high strength polyester yarn is double-knit in an interlocking pattern. The edges are sealed using a proprietary ultrasonic sealing process which improves partible and fiber retention, while remaining pyrogen-free.



APPLICATIONS

- These wipes are made from a continuous filament, double knit polyester which allows for consistent application of IPA to critical surfaces and equipment
- Great in areas where reduced IPA and volatile organic compounds may be necessary
- Their ultrasonic sealed edge improves fiber and particle retention for improved cleaning in critical environments
- Compatible with ISO Class 5 sterile cleanroom environments



SPECIFICATIONS

| Type | Method | Results |
|-----------|--------------------|---------------------------|
| Particles | Orbital Shake Test | 7×10^6 |
| Particles | IENT-4.3-6.1.4 | 0.2×10^6 |
| Fibers | IENT-4.3-6.2.2 | 230 fibers/m ² |
| Sterility | VD Max | >25 kGy |

| | |
|---------------|--|
| Part Number | W-BP7030PK1212S |
| Size | 12" X 12" |
| Packaging | 30 wipes/pack, 12 packs/case |
| Material | 100% pure polyester fiber |
| Solution | 70% IPA/30% WFI water |
| Basis Weight | 140 g/m ² |
| Compatibility | ISO Class 5 Sterile cleanroom environments |

V.7.BC

