

↑ Back to Table of Contents

- Connectors
- Ports & Flanges
- Bags
- Clamps & Clips
- Tubing
- Extension Lines
- Stopcocks
- Valves
- Injection & Sampling Ports

Filters

Guide Wire & Catheter Accessories

- Dilators, Introducer Sheaths & Needles
- Spikes
- Chambers
- Adhesives
- Syringes
- Sterilization Supplies
- Tools
- Applicators, Swabs & Brushes
- Drapes, Towels & Gauze
- Containers & Dispensers
- Clothing

FILTERS

Hydrophilic

Flow Direction →

40217
Hydrophilic Filter
NRFit™ Female Inlet
NRFit™ Male Outlet
Non-Vented Caps
0.2 Micron
PC, Nylon, PP

ISO
80369-6



ISO
80369-6

Also Available:	
Part No.	Material
40218	PC, PES, PP
40219	PC, PET, PP



12114
0.2 Micron
Both Ports 0.091 inch ID
(2.31 mm)
Acrylic, PES, PVDF



28216
Female Luer Lock Inlet
Male Luer Lock Outlet
0.2 Micron; 25 mm
Acrylic, Supor®



28217
Female Luer Lock Inlet
Male Luer Slip Outlet
0.2 Micron; 25 mm
Acrylic, Supor®



28245
Female Luer Lock Inlet
Male Luer Slip Outlet
0.2 Micron; 25 mm
PP, PES



28302
Minisart®
Female Luer Lock Inlet
Male Luer Slip Outlet
0.22 Micron; 28 mm
PES, MBS



28303
Minisart®
Female Luer Lock Inlet
Male Luer Lock Outlet
0.22 Micron; 28 mm
PES, MBS



28305
Female Luer Lock Inlet
Male Luer Lock Outlet
0.45 Micron; 25 mm
Acrylic, CA



28306
Female Luer Lock Inlet
Male Luer Lock Outlet
0.8 Micron; 25 mm
Acrylic, CA



28307
Female Luer Lock Inlet
Male Luer Lock Outlet
1.2 Micron; 25 mm
Acrylic, CA



28308
Female Luer Lock Inlet
Male Luer Lock Outlet
5 Micron; 25 mm
Acrylic, CA



28304
Female Luer Lock Inlet
Male Luer Lock Outlet
5 Micron; 25 mm
Acrylic, Supor®



28232
Filter Hub
Female Luer Lock Inlet
Male Luer Slip Outlet
5 Micron
ABS, Nylon



FILTERS

Hydrophilic

Flow Direction →



28243
Hydrophilic Filter Straw®
Female Luer Lock Outlet
5 Micron; 1.75 inch (44.45 mm) Straw
ABS, Nylon, DEHP-Free PVC



28244
Hydrophilic Filter Straw®
Female Luer Lock Outlet
5 Micron; 4 inch (101.6 mm) Straw
ABS, Nylon, DEHP-Free PVC



ISO
80369-6
40210
NRFit™ Filter Straw
5 Micron
1.75 inch (44.45 mm)
PVC, ABS, Nylon



ISO
80369-6
40211
NRFit™ Filter Straw
Non-Vented Cap
5 Micron
1.75 inch (44.45 mm)
PVC, ABS, Nylon, PP



ISO
80369-6
40212
NRFit™ Filter Straw
5 Micron
4 inch (101.6 mm)
PVC, ABS, Nylon



ISO
80369-6
40213
NRFit™ Filter Straw
Non-Vented Cap
5 Micron
4 inch (101.6 mm)
PVC, ABS, Nylon, PP

For more information about the
ISO 80369 series of standards,
including updates to ISO 80369-6,
please visit "Navigating ISO 80369" in
our Resource Center on qosina.com.

For more NRFit™ products, see
syringes on [pages 321](#) and
connectors and caps on [page 38](#).

Back to Table of Contents ↑

- Connectors
- Ports & Flanges
- Bags
- Clamps & Clips
- Tubing
- Extension Lines
- Stopcocks
- Valves
- Injection & Sampling Ports

Filters

Guide Wire & Catheter Accessories

- Dilators, Introducer Sheaths & Needles
- Spikes
- Chambers
- Adhesives
- Syringes
- Sterilization Supplies
- Tools
- Applicators, Swabs & Brushes
- Drapes, Towels & Gauze
- Containers & Dispensers
- Clothing