

POLYESTER WIPERS 9"X9"

SPECIFICATION SHEET



DESCRIPTION

The BenchmarkGMP Polyester Wipers are made from a 100% continuous-filament, double-knit polyester fabric. They are cleanroom packed and laundered with product traceability. Their ultrasonic sealed edge means they are extremely low particle and fiber generating. They are an ideal economic and quality wiper choice.



APPLICATIONS

- Ideal for ISO 4 cleanrooms and above
- Cleaning and disinfecting in all areas of cleanroom and production
- Cleanroom applications where abrasion and chemical resistance are required

SPECIFICATIONS

Type

Fibers
LPC
Organic Contaminates

Method

Orbital Shake Test IEST-4.3
Orbital Shake Test IEST-4.3

Results

< 400 counts/m²
6.0 x 10⁶/m²
Does not contain Silicon, Amides or DOP

| | |
|---------------------------------|--|
| Part Number | W-PKUW-99 |
| Size | 9"x 9" |
| Packaging | 150 wipes per bag, 10 bags per case |
| Material | 100% Virgin Grade Polyester (75 Denier) |
| Edge Model | Ultrasonic Sealed Edge |
| Basis Weight | 135 g/m ² |
| Sorptive Capacity | 380 ml/m ² |
| Sorptive Rate | < 1 second |
| Ions Contamination (ppm) | Calcium: 0.09 Sodium: 0.20 Potassium: 0.03 |
| NVR (g/m²) | DI-Water: ≤ 0.02 IPA: ≤ 0.03 |