

SINCE 1925

KEYSTONE

SAFETY ✦ CLEANROOM PRODUCTS

1422A TYVEK® EQUIPMENT COVERS



Tyvek® Cart and Rack Covers are designed to keep cleaned pharmaceutical equipment and supplies clean and protected.

TYPICAL APPLICATION & CHARACTERISTICS:

- Each cover is individually poly bagged and labeled for easy identification
- Available with, or without elastic; with or without zippers or snaps
- Breathable Tyvek® fabric allows residual moisture from the cleaning process to evaporate
- Hydrophobic
- Resists tearing and punctures
- Anti-static
- Lint-free
- Manufactured and packaged in ISO certified Cleanroom
- Certificates of Conformance sent with shipment
- Made in USA

MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES 1422A TYVEK®
Basis Weight	g/m2	EN ISO 536	41
Total Thickness	μ	EN 20534	130
Air Permeability	s	ISO 5636-5	20
Water Vapor Permeability	G/M2/24 HR	ASTM 398	1700
Hydrostatic Head	cm H2O	ISO 811	120
Trapezoidal Tear MD/CD	N	ISO 9073-4	26.1/30.6
Linting Propensity	Rating	BS 6909 (Shirley Method 21)	Pass Rating: Excellent
Burst Resistance	kPa	ISO 2960 (50 cm2)	108
Puncture Resistance	N	En 863	10.8
Melting Point	C°		135

SINCE 1925

KEYSTONE

SAFETY ✦ CLEANROOM PRODUCTS

1422A TYVEK® EQUIPMENT COVERS



KEYSTONE PART NUMBER	DIMENSIONS (L X W X H)
COV-TVK-141008E	14"X10"X8"
COV-TVK-262010E	26"X20"X10"
COV-TVK-301616E	30"X16"X16"
COV-TVK-302212E	30"X22"X12"
COV-TVK-302456E	30"X24"X56"
COV-TVK-322824E	32"X28"X24"
COV-TVK-362660E	36"X26"X60"
COV-TVK-423460E	42"X34"X60"
COV-TVK-442430E	44"X24"X30"
COV-TVK-443632E	44"X36"X32"
COV-TVK-444460E	44"X44"X60"
COV-TVK-482660E	48"X26"X60"
COV-TVK-482872E	48"X28"X72"
COV-TVK-484860E	48"X48"X60"
COV-TVK-502474E	50"X24"X74"
COV-TVK-502872E	50"X28"X72"
COV-TVK-505050E	50"X50"X50"
COV-TVK-545486E	54"X54"X86"
COV-TVK-563660E	56"X36"X60"
COV-TVK-602460E	60"X24"X60"
COV-TVK-602854E	60"X28"X54"
COV-TVK-603072E	60"X30"X72"
COV-TVK-604444E	60"X44"X44"
COV-TVK-605060E	60"X50"X60"
COV-TVK-642460E	64"X24"X60"
COV-TVK-642478E	64"X24"X78"
COV-TVK-643068E	64"X30"X68"
COV-TVK-644058E	64"X40"X58"
COV-TVK-682238E	68"X22"X38"
COV-TVK-742686E	74"X26"X86"
COV-TVK-762672E	76"X26"X72"