

SINCE 1925

KEYSTONE

SAFETY ✦ CLEANROOM PRODUCTS

1422A TYVEK® CYLINDER SHAPED EQUIPMENT COVERS



Cylinder shaped Tyvek® Covers are designed to keep cleaned pharmaceutical tanks, filters and reactors clean and protected. Also used for in-process protection to prevent cross-contamination of drug products.

TYPICAL APPLICATION & CHARACTERISTICS:

- Each cover is individually poly bagged and labeled for easy identification
- Available with, or without elastic
- Breathable Tyvek® fabric allows residual moisture from the cleaning process to evaporate
- Hydrophobic
- Resists tearing and punctures
- Anti-static
- Lint-Free
- Manufactured and packaged in ISO certified cleanrooms
- Certificates of Conformance sent with all shipments
- Made in USA

MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES 1422A TYVEK®
Basis Weight	g/m2	EN ISO 536	41
Total Thickness	μ	EN 20534	130
Air Permeability	s	ISO 5636-5	20
Water Vapor Permeability	G/M2/24 HR	ASTM 398	1700
Hydrostatic Head	cm H2O	ISO 811	120
Trapezoidal Tear MD/CD	N	ISO 9073-4	26.1/30.6
Linting Propensity	Rating	BS 6909 (Shirley Method 21)	Pass Rating: Excellent
Burst Resistance	kPa	ISO 2960 (50 cm2)	108
Puncture Resistance	N	En 863	10.8
Melting Point	C°		135

SINCE 1925

KEYSTONE

SAFETY ✦ CLEANROOM PRODUCTS

1422A TYVEK® CYLINDER SHAPED EQUIPMENT COVERS



KEYSTONE PART NUMBER	SIZE CHART
COV-TVK-1460E	14" Diameter X 60" Tall
COV-TVK-1613E	16" Diameter X 13" Tall
COV-TVK-2009E	20" Diameter X 9" Tall
COV-TVK-2048E	20" Diameter X 48" Tall
COV-TVK-3012E	30" Diameter X 12" Tall
COV-TVK-3260E	32" Diameter X 60" Tall
COV-TVK-3636E	36" Diameter X 36" Tall
COV-TVK-3648E	36" Diameter X 48" Tall
COV-TVK-4032E	40" Diameter X 32" Tall
COV-TVK-4410E	44" Diameter X 10" Tall
COV-TVK-4836E	48" Diameter X 36" Tall
COV-TVK-5050E	50" Diameter X 50" Tall
COV-TVK-5060E	50" Diameter X 60" Tall
COV-TVK-6060E	60" Diameter X 60" Tall
COV-TVK-7272E	72" Diameter X 72" Tall