

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

BLUE BHD FILLING NEEDLE BAG WITH TYVEK® DRAWSTRING



TYPICAL APPLICATION & CHARACTERISTICS:

- Designed for steam sterilization
- Incremental improvement compared to cutting and wrapping with Bio-Shield™ and autoclave tape
- Avoid sterilization Re-validation
- Tyvek® drawstring closure
- Steam sterilization indicator
- Eliminates need to cut material and the use of particulate generating autoclave tape
- Fast and consistent
- Creates secondary level of protection during aseptic presentation and setup of equipment
- Resists tearing and punctures
- 2-ply for good microbial barrier
- Manufactured and packaged in ISO certified Cleanrooms
- Certificates of Conformance sent with all shipments
- Made in USA

MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	BHD MATERIAL
Basis Weight	Gm/m2	57
Dry Tensile Strength (MD)	grams	8,100
Dry Tensile Strength (CD)	grams	7,100
Elongation Dry	%	27.3
Air Permeability	liters / 100 sq. cm	110
Suter Water Head (pressure test)	mm	400
Dry Mullen Burst	Grams/sq. cm	2,150
Flammability	seconds	8

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

BLUE BHD FILLING NEEDLE BAG WITH TYVEK® DRAWSTRING



KEYSTONE PART NUMBER	DIMENSIONS
BNWBG0104D-2P	1" X 4"
BNWBG0106D-2P	1" X 6"
BNWBG0108D-2P	1" X 8"
BNWBG0204D-2P	2" X 4"
BNWBG0206D-2P	2" X 6"
BNWBG0210D-2P	2" X 10"
BNWBG0308D-2P	3" X 8"
BNWBG0310D-2P	3" X 10"