

SINCE 1925

**KEYSTONE**

SAFETY ↔ CLEANROOM PRODUCTS

# BLUE BHD MUFFS WITH TYVEK® DRAWSTRING



Used as an autoclavable cover, these Blue Heavy Duty (BHD) Muffs eliminate the need for time consuming wrapping with sheets of Bioshield™ material and autoclave tape.

## TYPICAL APPLICATION & CHARACTERISTICS:

- Manufactured and packaged in ISO certified Cleanroom
- Drawstring closure allows for a tighter fit
- Autoclavable up to 135°C
- Conveniently packaged: 100 per heat-sealed, labeled poly bag
- Color: Blue
- Certificates of Conformance sent with all shipments
- Made in USA

## MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	BHD MATERIAL
Basis Weight	Gm/m2	57
Dry Tensile Strength (MD)	grams	8,100
Dry Tensile Strength (CD)	grams	7,100
Elongation Dry	%	27.3
Air Permeability	liters / 100 sq. cm	110
Suter Water Head (pressure test)	mm	400
Dry Mullen Burst	Grams / sq.cm	2,150
Flammability	seconds	8

SINCE 1925

**KEYSTONE**

SAFETY ↔ CLEANROOM PRODUCTS

# BLUE BHD MUFFS WITH TYVEK® DRAWSTRING



KEYSTONE PART NUMBER	DIMENSIONS
MFNWB05D	5" used for openings 1" to 2 1/2"
MFNWB07D	7" used for openings 1 1/2" to 4"
MFNWB08D	8" used for openings 2 1/2" to 4 3/4"
MFNBW09D	9" used for openings 3" to 5 1/2"
MFNWB12D	12" used for openings 4 1/2" to 8"
MFNWB16D	16" used for openings 7" to 12"
MFNWB18D	18" used for openings 9" to 14"