

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

HEAVY DUTY SMS AUTOCLAVE BAG, SELF-SEAL, WITH STEAM INDICATOR



TYPICAL APPLICATION & CHARACTERISTICS:

- Designed for higher temperature steam sterilization cycles (135°C)*
- Significant improvement versus cutting and wrapping with sheets and tape
- Avoid sterilization revalidation if currently wrapping with SMS sheets
- Manufactured and packaged in a certified ISO6 Cleanroom
- Clean--no cellulose
- Heat Sealed seams
- Convenient, time saving self-seal closure system
- Strong, rip and tear resistant 50 gsm SMS material
- Compatible with ethylene oxide and hydrogen peroxide sterilization
- An extensive range of sizes
- Excellent barrier characteristics prevent liquid and microbial penetration
- Sub packed into sealed labeled poly bags with manufacturing lot number
- Date of Manufacture on label for complete traceability
- Certificate of Conformance sent with each shipment
- Made in the USA

*Tyvek® is not suitable for steam sterilization cycles above 127°C.

MATERIAL SPECIFICATIONS:

PROPERTY	SPECIFICATION	UNIT
Basis Weight	50	Gm/m2
Grab Strength (CD)	20	lbs
Tear Strength (MD)	7.5	lbs
Mullen Burst Strength (Gross)	34	psi
Air Permeability	50	cfm
Hydrohead	55	mbar

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

HEAVY DUTY SMS AUTOCLAVE BAG, SELF-SEAL, WITH STEAM INDICATOR



KEYSTONE PART NUMBER	DIMENSIONS	KEYSTONE PART NUMBER	DIMENSIONS
BGSMS3SS09.516	9.5" x 16" 24 cm x 41 cm	BGSMS3SS2126	21" x 26" 53 cm x 66 cm
BGSMS3SS1217	12" x 17" 30 cm x 43 cm	BGSMS3SS2424	24" x 24" 61 cm x 61 cm
BGSMS3SS1221	12" x 21" 30 cm x 53 cm	BGSMS3SS2430	24" x 30" 61 cm x 76 cm
BGSMS3SS1318	13" x 18" 33 cm x 46 cm	BGSMS3SS2435	24" x 35" 61 cm x 89 cm
BGSMS3SS14.522	14.5" x 22" 37 cm x 56 cm	BGSMS3SS2630	26" x 30" 66 cm x 76 cm
BGSMS3SS1520	15" x 20" 38 cm x 51 cm	BGSMS3SS2836	28" x 36" 71 cm x 91 cm
BGSMS3SS1622	16" x 22" 41 cm x 56 cm	BGSMS3SS3030	30" x 30" 76 cm x 76 cm
BGSMS3SS1824	18" x 24" 46 cm x 61 cm	BGSMS3SS3636	36" x 36" 91 cm x 91 cm
BGSMS3SS1830	18" x 30" 46 cm x 76 cm	BGSMS3SS3648	36" x 48" 91 cm x 122 cm
BGSMS3SS2028	20" x 28" 51 cm x 71 cm	BGSMS3SS4848	48" x 48" 122 cm x 122 cm