

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

SMS STERILIZATION WRAPPER HEAVY DUTY 50 GSM



Keystone Cleanroom Products' line of Heavy Duty SMS wrappers provide an excellent microbial barrier. Ideal for pharmaceutical production applications.

- Manufactured and packaged in a certified ISO8 Cleanroom
- Clean--no cellulose
- Strong, rip and tear resistant 50 gsm SMS material
- Soft and drape-able fabric
- An extensive range of sizes available means less cutting
- Hydrophobic barrier allows sterilization gasses, but prevents liquid and microbial penetration
- Sub packed into sealed labeled poly bags with Lot Number and Date of Manufacture for complete traceability
- Certificate of Conformance sent with each shipment
- Made in the USA

For use in steam, ethylene oxide and hydrogen peroxide sterilization processes.

MATERIAL SPECIFICATIONS:


| PROPERTY | SPECIFICATION | UNIT |
|-------------------------------|---------------|-------------------|
| Basis Weight | 50 | Gm/m ² |
| Grab Strength (CD) | 20 | lbs |
| Tear Strength (MD) | 7.5 | lbs |
| Mullen Burst Strength (Gross) | 34 | psi |
| Air Permeability | 50 | cfm |
| Hydrohead | 55 | mbar |

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

SMS STERILIZATION WRAPPER HEAVY DUTY 50 GSM



| KEYSTONE PART NUMBER | DIMENSIONS | KEYSTONE PART NUMBER | DIMENSIONS |
|----------------------|-------------|----------------------|------------|
| KCRP-300-01.5X01.5 | 1.5" X 1.5" | KCRP-300-2020 | 20" X 20" |
| KCRP-300-0202 | 2" X 2" | KCRP-300-2424 | 24" X 24" |
| KCRP-300-0303 | 3" X 3" | KCRP-300-3030 | 30" X 30" |
| KCRP-300-0404 | 4" X 4" | KCRP-300-3636 | 36" X 36" |
| KCRP-300-0606 | 6" X 6" | KCRP-300-4040 | 40" X 40" |
| KCRP-300-0808 | 8" X 8" | KCRP-300-4545 | 45" X 45" |
| KCRP-300-0909 | 9" X 9" | KCRP-300-4548 | 45" X 48" |
| KCRP-300-0912 | 9" X 12" | KCRP-300-4848 | 48" X 48" |
| KCRP-300-1212 | 12" X 12" | KCRP-300-5454 | 54" X 54" |
| KCRP-300-1515 | 15" X 15" | KCRP-300-5472 | 54" X 72" |
| KCRP-300-1818 | 18" X 18" | KCRP-300-6060 | 60" X 60" |
| KCRP-300-1822 | 18" X 22" | KCRP-300-6080 | 60" X 80" |