

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

BLUE SUPER DUTY SMS STOPPER BOWL COVER



Save time by using the stopper bowl cover instead of wrapping with sheets of SMS material and autoclave tape. When used in conjunction with a secondary autoclave bag, it creates an effective wrapping system. After the outer bag is removed, this stopper bowl cover should remain in place, providing protection from microbial and particulate contamination while the remaining set-up activities are performed.

TYPICAL APPLICATION & CHARACTERISTICS:

- Eliminates autoclave tape residue
- SMS has excellent microbial barrier properties
- Can be sterilized in **135°** autoclave cycles
- Creates secondary level of protection during aseptic presentation and setup of equipment
- Hydrophobic
- Resists tearing and punctures
- Lint-free
- Manufactured and packaged in ISO certified cleanrooms
- Conveniently packaged: 1 per heat-sealed and labeled poly bag
- Certificates of Conformance sent with all shipments
- Made in USA

MATERIAL SPECIFICATIONS:


PROPERTY	SPECIFICATION	UNIT
Basis Weight	88	Gm/m ²
Grab Strength (CD)	40	lbs
Tear Strength (MD)	11	lbs
Mullen Burst Strength (Gross)	60	psi
Air Permeability	25	cfm
Hydrohead	70	cm

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

BLUE SUPER DUTY SMS STOPPER BOWL COVER



KEYSTONE PART NUMBER	DIMENSIONS
BWL-KCRP-600-20	Fits 11" to 16" Diameter
BWL-KCRP-600-24	Fits 12" to 20" Diameter
BWL-KCRP-600-26	Fits 14" to 22" Diameter
BWL-KCRP-600-30	Fits 20" to 26" Diameter
BWL-KCRP-600-36	Fits 26" to 30" Diameter
BWL-KCRP-600-42	Fits 28" to 35" Diameter
BWL-KCRP-600-46	Fits 30" to 38" Diameter
BWL-KCRP-600-48	Fits 30" to 40" Diameter
BWL-KCRP-600-54	Fits 35" to 45" Diameter
BWL-KCRP-600-60	Fits 40" to 50" Diameter