

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

UV LIGHT BLOCKING TYVEK® FC EQUIPMENT COVERS WITH ZIPPERS



UV Light Blocking Tyvek® FC Equipment Covers are designed to keep cleaned pharmaceutical equipment and supplies clean and protected from external sources of contamination while also providing in-process protection to light sensitive biotech products.

TYPICAL APPLICATION & CHARACTERISTICS:

- Each cover is individually poly bagged and labeled for easy identification
- Available with, or without elastic with, or without zippers, or snaps
- UV Light Blocking Tyvek® FC fabric protects light sensitive Biotech products
- Hydrophobic
- Resists tearing and punctures
- Anti-static
- Manufactured and packaged in ISO certified Cleanroom
- Certificates of Conformance sent with each shipment
- Made in USA

MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES TYVEK® 400 FC
Basis Weight	oz/yd ²	ASTM D3776	1.6 oz/yd ²
Thickness	mils	ASTM D1117	7.4 mils
Breaking Strength-Grab (MD)	f	ASTM D5034	20 lbf
Breaking Strength-Grab (CD)	f	ASTM D5034	26 lbf
Hydrostatic Head	in H ₂ O	AATCC 127	4.2 in H ₂ O
Tear Resistance MD	f	ASTM D5733	5lbf
Tear Resistance CD	f	ASTM D5733	7lbf
Flammability		16 CFR 1610	Class 1

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

UV LIGHT BLOCKING TYVEK® FC COVERS WITH ZIPPERS



KEYSTONE PART NUMBER	DIMENSIONS
363235TGZ4	36" x 32" x 35"
363268TGZ4	36" x 32" x 68"
563666TGZ2	56" x 36" x 66"
X604456TGZ2WP2	60" x 44" x 56"
604860Z2TGWP2	60" x 48" x 60"